

HOLE-TYPE OR LAMP OPL39-50



HOLE-TYPE OR LAMP OPL39-50

Hole type operation light is suitable for small operating room. Labstac lights uses multiple different sources and by focusing knob, the beam is illuminated from different surgery locations which reduces the shadow by medical staff head, hands and other parts caused. The illumination adjusted according to actual needs.

Used in Medical unit, Surgery, Operation theatre.

Also known as Hole type surgical Light, Hole type OT light.

OPL39-50 HOLE-TYPE OR LAMP



SPECIFICATIONS

Model	OPL39-50
Illumination	40000 Lux
The value of the total irradiance	$<6\text{mW} / \text{m}^2\text{lx}$
Total irradiance	$\leq 600\text{W}/\text{m}^2$
Color rendering Index (Ra)	$100 \geq \text{Ra} \geq 85$
Color Temperature (K)	$4000 \pm 500 \text{ K}$
Power supply	$220 \pm 10\% \text{ 50 Hz DC24V}$
Power of bulb	100 W
Input power	180 VA



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com