

# HOLE-TYPE OR LAMP OPL39-50I



# HOLE-TYPE OR LAMP OPL39-50I

Hole type operation light is suitable for small operating room. Labstac lights uses multiple different sources and by focusing knob, the beam is illuminated from different surgery locations which reduces the shadow by medical staff head, hands and other parts caused. The illumination adjusted according to actual needs.

Used in Medical unit, Surgery, Operation theatre.

Also known as Hole type surgical Light, Hole type OT light.

## OPL39-50I HOLE-TYPE OR LAMP



### SPECIFICATIONS

Model	OPL39-50I
Diameter of lighthouse (cm)	Combined type
Illumination at 1m	40000-100000 Lux
Light field diameter at 1m (cm)	10-28
D10 (mm)	180±40
Diameter of spot D50 (mm)	110±35
L1+L2 Depth of illumination	150 cm
Color rendering Index (Ra)	100≥Ra≥85
Color Temperature (K)	4500±500
Number of LED's	32
Average bulb life (hours)	50000 h
Heat-to-light ratio (mW/m <sup>2</sup> lx)	<6
Working distance	70-140 cm
Input power	60 VA
Ceiling height for installation (cm)	2.7-3.2 m



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)