

# HOLE-TYPE OR LAMP OPL3B-12



# HOLE-TYPE OR LAMP OPL3B-12

Hole type operation light is suitable for small operating room. Labstac lights uses multiple different sources and by focusing knob, the beam is illuminated from different surgery locations which reduces the shadow by medical staff head, hands and other parts caused. The illumination adjusted according to actual needs.

Used in Medical unit, Surgery, Operation theatre.

Also known as Hole type surgical Light, Hole type OT light.

## OPL3B-12 HOLE-TYPE OR LAMP



### SPECIFICATIONS

|                                   |                                      |
|-----------------------------------|--------------------------------------|
| Model                             | OPL3B-12                             |
| Illumination                      | 50000 Lux                            |
| Depth of light beam               | $\geq 500$ mm                        |
| The value of the total irradiance | $< 6\text{mW} / \text{m}^2\text{lx}$ |
| Total irradiance                  | $\leq 600\text{W}/\text{m}^2$        |
| Color rendering Index (Ra)        | $100 \geq \text{Ra} \geq 85$         |
| Color Temperature (K)             | $4000 \pm 500$ K                     |
| Power supply                      | $220 \pm 10\%$ 50 Hz DC24V           |
| Power of bulb                     | 125 W                                |
| Input power                       | 150 VA                               |



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)