

# MANUAL OR TABLE OPT12-2000



# MANUAL OR TABLE OPT12-2000

Manual Operation table are suitable for various operation tables conducted by medical unit. The OR table and accessories are ergonomic and can be used across multiple disciplines. Various actions are controlled by hydraulic pneumatic system. 360° rotation helps every surgical requirement.

Used in Medical unit, Surgery, Operation theatre, cardiovascular, pediatric, gynecology, abdominal, gall bladder, plastic surgery and more.

Also known as Manual surgical table, Manual OR Table, Manual operating table.

## OPT12-2000 MANUAL OR TABLE

Imported LINAK motor various actions are controlled by remote controller

The table, base seat and upright column covering are made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned

The mattress adopts high density memory sponge and is once formed, with no seam

Whole stainless steel hidden auxiliary table with flexible extension and retraction

Equipped with stainless steel sewage basin and other accessories

Foot treadle brake device which is easy for moving obstetric table



## SPECIFICATIONS

Model	OPT12-2000
Table Length	2000±50 mm
Table Width	600±20 mm
Height Adjustment	(700~950) ±50 mm
Working area length	520 mm
Working area width	500 mm
Trendelenburg	≥10°
Reverse trendelenburg	≥12°
Backboard upwards	≥50°



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)