

MANUAL OR TABLE OPT13-2000



MANUAL OR TABLE OPT13-2000

Manual Operation table are suitable for various operation tables conducted by medical unit. The OR table and accessories are ergonomic and can be used across multiple disciplines. Various actions are controlled by hydraulic pneumatic system. 360° rotation helps every surgical requirement.

Used in Medical unit, Surgery, Operation theatre, cardiovascular, pediatric, gynecology, abdominal, gall bladder, plastic surgery and more.

Also known as Manual surgical table, Manual OR Table, Manual operating table.

OPT13-2000 MANUAL OR TABLE

Various actions are controlled by hydraulic pneumatic system and mechanical driven system

Head plate can be exchanged with leg plate which can fulfill the needs of C type arm complete examination, X-ray photography

Imported Y type sealing ring is durable

The table surface adopts quality carbon plastic plate, the mattress adopts high density memory sponge and is once formed, with no seam

Quality stainless steel 304 is used for base plate, column shield and rails on both sides of the table, easy for cleaning and sterilization

Foot treadle brake device which is easily for moving operating table



SPECIFICATIONS

Model	OPT13-2000
Table Length	2000±50 mm
Table Width	500±20 mm
Height Adjustment	(700~950) ±50 mm
Left tilt	≥15°
Right tilt	≥15°
Trendelenburg	≥15°
Reverse trendelenburg	≥15°
Head board upwards	≥30°
Head board downwards	≥90°
Backboard upwards	≥60°
Backboard downwards	≥20°
Leg board downwards	≥90°



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com