

DRYING OVEN OVE11-070



DRYING OVEN OVE11-070

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

OVE11-070 DRYING OVEN

Precision LCD

Cooling fan

Stainless steel inner chamber

Intelligent PID control system

Best quality oven



SPECIFICATIONS

Model	OVE11-070
Capacity	70 L
Temperature Range	RT+5-250°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±1%
Temperature Uniformity	±1°C
Temperature Resolution	±0.5°C
Timer Range	0-999 min
Internal Dimension (WxDxH)	450x350x420 mm
Overall Dimension (WxDxH)	740x550x630 mm
Package Size	720x740x1470 mm
Weight	47/50 kg
Power	1200 W
Power Supply	220 V, 50 Hz



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com