

# DRYING OVEN OVE11-070



# DRYING OVEN OVE11-070

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

## OVE11-070 DRYING OVEN

Precision LCD

Cooling fan

Stainless steel inner chamber

Intelligent PID control system

Best quality oven



## SPECIFICATIONS

|                            |                 |
|----------------------------|-----------------|
| Model                      | OVE11-070       |
| Capacity                   | 70 L            |
| Temperature Range          | RT+5-250°C      |
| Temperature Accuracy       | ±0.1°C          |
| Temperature Fluctuation    | ±1%             |
| Temperature Uniformity     | ±1°C            |
| Temperature Resolution     | ±0.5°C          |
| Timer Range                | 0-999 min       |
| Internal Dimension (WxDxH) | 450x350x420 mm  |
| Overall Dimension (WxDxH)  | 740x550x630 mm  |
| Package Size               | 720x740x1470 mm |
| Weight                     | 47/50 kg        |
| Power                      | 1200 W          |
| Power Supply               | 220 V, 50 Hz    |



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)