

# DRYING OVEN OVE11-225



# DRYING OVEN OVE11-225

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

## OVE11-225 DRYING OVEN

Precision LCD

Cooling fan

Stainless steel inner chamber

Intelligent PID control system

Best quality oven



## SPECIFICATIONS

Model	OVE11-225
Capacity	225 L
Temperature Range	RT+5-250°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±1%
Temperature Uniformity	±1°C
Temperature Resolution	±0.5°C
Timer Range	0-999 min
Convection Mode	Forced Convection
Internal Dimension (WxDxH)	600x550x750 mm
Overall Dimension (WxDxH)	965x770x980 mm
Package Size	1907x815x1090 mm
Shelves	4
Weight	96/105 kg
Power	2500 W
Power Supply	220 V, 50 Hz



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)