

# DRYING OVEN OVE12-040



# DRYING OVEN OVE12-040

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

## OVE12-040 DRYING OVEN

Digital display all 304 stainless steel inner chamber and outer

Intelligent PID control system

Large observation window

Tight buckle handle



## SPECIFICATIONS

Model	OVE12-040
Capacity	40 L
Temperature Accuracy	$\pm 1^{\circ}\text{C}$
Temperature Fluctuation	$\pm 1\%$
Temperature Uniformity	$0.5^{\circ}\text{C}$
Temperature Resolution	$0.1^{\circ}\text{C}$
Time for Heating	$\leq 6$ min ( $20^{\circ}\text{C}$ to $57^{\circ}\text{C}$ )
Timer Range	0-999 min
Internal Dimension (WxDxH)	350×350×350 mm
Overall Dimension (WxDxH)	640×550×530 mm
Weight	35.5/38.5 kg
Power	1000 W
Power Supply	220 V, 50 Hz



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)