

DRYING OVEN OVE12-070



DRYING OVEN OVE12-070

Heating Oven, Laboratory Drying Oven

Used in Drying, Roasting, Laboratory, Research, Engineering, Industry, Agricultural, Biology.

Also known as Laboratory Oven, Laboratory general Purpose Oven.

OVE12-070 DRYING OVEN

Micro-computer intelligent control system with LED display

Forced air convection maintains temperature uniformity

High quality stainless steel inner chamber with electrostatic spray treated exterior

Large observation window

Easy to clean semicircular corners with adjustable shelves

Silicon rubber sealing ring prevents heat loss

Over temperature alarm with timer



SPECIFICATIONS

Model	OVE12-070
Capacity	70 L
Temperature Range	RT+5-250°C
Temperature Accuracy	±1°C
Temperature Uniformity	±4.5°C at 250°C
Temperature Resolution	±1°C
Timer Range	0-999 min
Internal Dimension (WxDxH)	420x450x350 mm
Overall Dimension (WxDxH)	630x740x550 mm
Package Size	750x820x630 mm
Shelves	2
Weight	44.5/47.5 kg
Power	1200 W
Power Supply	220V, 60Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS54389	Intelligent program temperature controller	
LS54401	Independent limit of temperature controller	



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com