

# THERMOSTATIC DRYING OVEN OVE13-040



# THERMOSTATIC DRYING OVEN OVE13-040

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

## OVE13-040 THERMOSTATIC DRYING OVEN

ALLSENSTM programmable PID control

ALLHEATTM cavity warm-up technology

ALLFLOWTM perfect air current cycling

Pro-INSULATIONTM Insulation technology

SUS304 stainless steel internal material

Build-in memory menu, real-time monitoring

Display parameters: Temperature, Air velocity, Run time, Run/Stop



## SPECIFICATIONS

| Model                      | OVE13-040                |
|----------------------------|--------------------------|
| Capacity                   | 40 L                     |
| Temperature Range          | RT+5 – 250°C             |
| Ambient Temperature        | 10-30°C                  |
| Temperature Accuracy       | 0.1°C                    |
| Temperature Fluctuation    | ±1°C(50 - 240°C)         |
| Temperature Uniformity     | 2%(50 - 240°C)           |
| Ambient Humidity           | <70%                     |
| Timer Range                | 0-99 h or 0-9999 min     |
| Convection Mode            | Forced Convection        |
| Internal Dimension (WxDxH) | 350x350x350 mm           |
| Overall Dimension (WxDxH)  | 570x580x593 mm           |
| Package Size               | 690x660x680 mm           |
| Shelves                    | 2                        |
| Weight                     | 40/43 kg                 |
| Power                      | 770 W                    |
| Power Supply               | Single phase AC220V/50Hz |





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)