

# VACUUM DRYING OVEN OVE21-25S



# VACUUM DRYING OVEN OVE21-25S

Vacuum ovens are very versatile pieces of equipment with applications in laboratory research, engineering, and industry. The low pressure environment also minimizes oxidation during drying. A rugged, high quality pressure chamber, gasket seals, and pump are features to look for, as well as convenient, programmable controls and a computer interface.

Used in Biochemical, Pharmaceutical, Healthcare, Agricultural, Environmental protection, disinfection, Sterilization, Drying, Roasting, Laboratory, Research, engineering, Industry,.

Also known as Vacuum Oven, Laboratory Vacuum Oven.

## OVE21-25S VACUUM DRYING OVEN

- Stainless steel inner chamber without scroll rolled
- Silicone rubber door seal
- Strong sealing
- Nice design



### SPECIFICATIONS

Model	OVE21-25S
Capacity	25 L
Temperature Accuracy	±1°C
Temperature Fluctuation	±2%
Vacuum Degree	< 133 Pa
Timer Range	0-999 min
Internal Dimension (WxDxH)	300x300x280 mm
Shelves/Trays	2
Power	900 W
Power Supply	220 V, 50 Hz



#### Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)