

VACUUM DRYING OVEN OVE21-25S



VACUUM DRYING OVEN OVE21-25S

Vacuum ovens are very versatile pieces of equipment with applications in laboratory research, engineering, and industry. The low pressure environment also minimizes oxidation during drying. A rugged, high quality pressure chamber, gasket seals, and pump are features to look for, as well as convenient, programmable controls and a computer interface.

Used in Biochemical, Pharmaceutical, Healthcare, Agricultural, Environmental protection, disinfection, Sterilization, Drying, Roasting, Laboratory, Research, engineering, Industry,.

Also known as Vacuum Oven, Laboratory Vacuum Oven.

OVE21-25S VACUUM DRYING OVEN

Stainless steel inner chamber without scroll rolled

Silicone rubber door seal

Strong sealing

Nice design



SPECIFICATIONS

Model	OVE21-25S
Capacity	25 L
Temperature Accuracy	$\pm 1^{\circ}\text{C}$
Temperature Fluctuation	$\pm 2\%$
Vacuum Degree	$< 133 \text{ Pa}$
Timer Range	0-999 min
Internal Dimention (WxDxH)	300x300x280 mm
Shelves/Trays	2
Power	900 W
Power Supply	220 V, 50 Hz



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com