

PIPETTE TIPS

PIC11-05000



PIPETTE TIPS PIC11-05000

Pipette tips are essential when using a micropipette. These get attached to the shaft of the pipette and are usually ejected after use. It is an essential lab item for any test or experiment. They are available as racked or loose bulk tips, with or without a filter, sterile or non-sterile.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Pipet Tips, Laboratory Pipette Tips, Laboratory Pipet Tips.

PIC11-05000 PIPETTE TIPS

Made from medical grade Polypropylene

Slippy surface with excellent design

DNase-free

Rnase-free

Non-pyrogenic

Filter and non-filter for choice

Bag and rack are available to pack



SPECIFICATIONS

| | |
|---------|--------------|
| Model | PIC11-05000 |
| Volume | 5000 μ l |
| Package | 100 pcs/bag |
| Unit | 10 |



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com