

ULTRA LOW RETENTION TIP PIC12-01000



ULTRA LOW RETENTION TIP PIC12-01000

Pipette tips are essential when using a micropipette. These get attached to the shaft of the pipette and are usually ejected after use. It is an essential lab item for any test or experiment. They are available as racked or loose bulk tips, with or without a filter, sterile or non-sterile.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Pipet Tips, Laboratory Pipette Tips, Laboratory Pipet Tips.

PIC12-01000 ULTRA LOW RETENTION TIP

Made from medical grade Polypropylene

Slippy surface with excellent design


DNase-free

Rnase-free

Non-pyrogenic

Filter and non-filter for choice

Bag and rack are available to pack



SPECIFICATIONS

| | |
|---------|--------------|
| Model | PIC12-01000 |
| Volume | 1000 μ l |
| Package | 1000 pcs/bag |
| Unit | 10 |



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com