

# RACKED ULTRA LOW RETENTION TIP 96X5 PIC17-0010T1



# RACKED ULTRA LOW RETENTION TIP 96X5

## PIC17-0010T1

Pipette tips are essential when using a micropipette. These get attached to the shaft of the pipette and are usually ejected after use. It is an essential lab item for any test or experiment. They are available as racked or loose bulk tips, with or without a filter, sterile or non-sterile.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Pipet Tips, Laboratory Pipette Tips, Laboratory Pipet Tips.

## PIC17-0010T1 RACKED ULTRA LOW RETENTION TIP 96X5

Made from medical grade Polypropylene

Slippy surface with excellent design

DNase-free

Rnase-free

Non-pyrogenic

Filter and non-filter for choice

Bag and rack are available to pack



### SPECIFICATIONS

Model	PIC17-0010T1
Volume	10 $\mu$ l
Unit	10



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)