

# RACKED ULTRA LOW RETENTION TIP 96X10 PIC17-1000T2



# RACKED ULTRA LOW RETENTION TIP 96X10

## PIC17-1000T2

Pipette tips are essential when using a micropipette. These get attached to the shaft of the pipette and are usually ejected after use. It is an essential lab item for any test or experiment. They are available as racked or loose bulk tips, with or without a filter, sterile or non-sterile.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Pipet Tips, Laboratory Pipette Tips, Laboratory Pipet Tips.

## PIC17-1000T2 RACKED ULTRA LOW RETENTION TIP 96X10

Made from medical grade Polypropylene

Slippy surface with excellent design

DNase-free

Rnase-free

Non-pyrogenic

Filter and non-filter for choice

Bag and rack are available to pack



### SPECIFICATIONS

Model	PIC17-1000T2
Volume	1000 $\mu$ l
Unit	10



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)