

# PIPETTOR SOLUTION BASINS PIC21-50



# PIPETTOR SOLUTION BASINS PIC21-50

Basins are designed for single and multi channel pipettors. Easy access to solution. It avoids spills and provides stability. They are available pre-sterile or nonsterile and constructed of either polyvinyl chloride or polystyrene. Commonly used in laboratories, research, pharmaceutical, etc.

Used in Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research.

Also known as Laboratory Pipettor Solution Basins.

## PIC21-50 PIPETTOR SOLUTION BASINS

Permits easy access to solution

Easy Cleaning



### SPECIFICATIONS

Model	PIC21-50
Volume	50 ml



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)