

# STACKABLE 1000ML TIP BOX PIC67-01000S



# STACKABLE 1000ML TIP BOX PIC67-01000S

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology. Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

## PIC67-01000S STACKABLE 1000ML TIP BOX

Made from medical grade Polypropylene  
Slippy surface with excellent design  
DNase-free  
Rnase-free  
Non-pyrogenic  
Filter and non-filter for choice  
Bag and rack are available to pack



## SPECIFICATIONS

Model	PIC67-01000S
Volume	1000 $\mu$ l
Unit	20



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)