

# SEROLOGICAL PIPETTE PIC71-3200C2



# SEROLOGICAL PIPETTE PIC71-3200C2

Manufactured for accurate and rapid fluid dispensing. It is sterile, pyrogen-free, cotton plugged, and can be used for a variety of liquid handling applications. The reservoir eliminates the need for multiple reagent containers and saves bench top space.

Used in Plasmid Purification, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Laboratory Serological Pipette.

## PIC71-3200C2 SEROLOGICAL PIPETTE

Easy to use

Medical grade original Polystyrene


RNase- free, DNase- free and Non-pyrogenic

Seven types of volume 1ml, 2ml, 5ml, 10ml, 25ml, 50ml and 100ml

Bidirectional graduation improve accuracy

Seven types of color code on tips make pipette

Sterilized by Gamma irradiation



## SPECIFICATIONS

Model	PIC71-3200C2
Volume	10 ml
Precision	2 %
Negative Graduation (ml)	-3
Graduation	1 / 10
Color	Red
Sterile	Gamma
Package	Bulk pack
Pcs/bag	25 pcs/bag
Pcs/ctn	3200 pcs/carton
Unit	1



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)