

SEROLOGICAL PIPETTE PIC71-4000C1



SEROLOGICAL PIPETTE PIC71-4000C1

Manufactured for accurate and rapid fluid dispensing. It is sterile, pyrogen-free, cotton plugged, and can be used for a variety of liquid handling applications. The reservoir eliminates the need for multiple reagent containers and saves bench top space.

Used in Plasmid Purification, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Laboratory Serological Pipette.

PIC71-4000C1 SEROLOGICAL PIPETTE

Easy to use

Medical grade original Polystyrene


RNase- free, DNase- free and Non-pyrogenic

Seven types of volume 1ml, 2ml, 5ml, 10ml, 25ml, 50ml and 100ml

Bidirectional graduation improve accuracy

Seven types of color code on tips make pipette

Sterilized by Gamma irradiation



SPECIFICATIONS

Model	PIC71-4000C1
Volume	1 ml
Precision	2 %
Negative Graduation (ml)	-0.4
Graduation	1 / 100
Color	Yellow
Sterile	Gamma
Package	Individual
Pcs/bag	100 pcs/bag
Pcs/ctn	4000 pcs/carton
Unit	1



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com