

# VARIABLE VOLUME 12 CHANNEL MICROPIPETTE PIP22-300E



# VARIABLE VOLUME 12 CHANNEL MICROPIPETTE

## PIP22-300E

New ergonomic technology and innovative design makes it one of the most innovative pipettes in the market. Multi Channel Pipettes are designed to aspirate and dispense liquid volumes while eliminating many of the potential error risks. Its operation remains correct and controlled, volume settings can be set to prevent accidental volume changes. Used in ELISA, PCR, Cell Culture, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology.

Also known as Laboratory Multi Channel Pipette.

## PIP22-300E VARIABLE VOLUME 12 CHANNEL MICROPIPETTE

Light weight

Rotating dispensing head for optimal pipetting convenience.

Due to the optimized standard design of the tip cone, it can be compatible with a wide range of hand of tips.

Large central pipetting button and separate ejection function

True one-handed operation for both right-and left-handers

Completely autoclavable at 121°C(20min)

Volume-change protection

4-Position volume display, always clearly visible externally

Short stroke of only 12.5mm to reduce the risk of RSI( Repetitive Strain Injury)

Corrosion-resistant piston and ejector

Color-coded for easy selection of the right tip

Individual shafts with seals can be easily unscrewed for cleaning or replacing-eliminates expense and long outages.



## SPECIFICATIONS

Model	PIP22-300E
Autoclavable	Ejector 121°C
Volume	30-300 µl
Volume Step (µl)	300/150/30
Increment (µl)	0.3
Accuracy* $\leq \pm\%$	0.6/1.2/3
CV* $\leq\%$	0.3/0.6/1.5
Type of tips (µl)	300
Unit	5



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)