

ADJUSTABLE GLASS- INJECTION DISPENSER AMBER GLASS PIP51-05AG



ADJUSTABLE GLASS-INJECTION DISPENSER

AMBER GLASS PIP51-05AG

A non-fragile, silicone tube manufactured from special materials. It is conveniently used with high efficiency. It is fully autoclavable and mostly used in laboratory, pharmaceutical and for other research purpose.

Used in Chemical Handling, Pollution Testing, Soil Testing, Environmental Laboratories, Pharmaceutical Laboratories, Research and Development, Qualitative and Quantitative Analysis.

Also known as Laboratory Liquid Dispenser.

PIP51-05AG ADJUSTABLE GLASS-INJECTION DISPENSER AMBER GLASS

Availability for autoclave of the whole device

A non-fragile, durable silicone tube made of special materials

Convenience of use with high efficiency, and low repetitive

Error, capable of greatly reducing time and costs in research and experiment

Solvent Resistance: general acid and alkali

Material: PP



SPECIFICATIONS

Model	PIP51-05AG
Volume	0.5-5 ml
Test Volume	0.5 ml
Accuracy	±0.04
Gross Weight (kg)	4.2
Dimensions in cubic meters	0.62x0.25x0.315
Package	10 pcs/carton



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com