

ADJUSTABLE GLASS- INJECTION DISPENSER TRANSLUCENT GLASS PIP51-10TG



ADJUSTABLE GLASS-INJECTION DISPENSER

TRANSLUCENT GLASS PIP51-10TG

A non-fragile, silicone tube manufactured from special materials. It is conveniently used with high efficiency. It is fully autoclavable and mostly used in laboratory, pharmaceutical and for other research purpose.

Used in Chemical Handling, Pollution Testing, Soil Testing, Environmental Laboratories, Pharmaceutical Laboratories, Research and Development, Qualitative and Quantitative Analysis.

Also known as Laboratory Liquid Dispenser.

PIP51-10TG ADJUSTABLE GLASS-INJECTION DISPENSER

TRANSLUCENT GLASS

Availability for autoclave of the whole device

A non-fragile, durable silicone tube made of special materials

Convenience of use with high efficiency, and low repetitive

Error, capable of greatly reducing time and costs in research and experiment

Solvent Resistance: general acid and alkali

Material: PP



SPECIFICATIONS

Model	PIP51-10TG
Volume	1-10 ml
Test Volume	1 ml
Accuracy	±0.1
Gross Weight (kg)	6.5
Dimensions in cubic meters	0.66x0.315x0.335
Package	10 pcs/carton



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com