

# 2°C-10°C MEDICAL REFRIGERATOR REF11-300



# 2°C-10°C MEDICAL REFRIGERATOR REF11-300

Engineered to meet the demanding requirements of laboratory research. Major feature like forced air circulation system makes this product highly effective and reliable to provide quick freezing. Designed Environment friendly with no ozone-damaging chemicals and space saving with easier cleanability.

Used in Industry, Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

Also known as Laboratory Medical Refrigerator.

## REF11-300 2°C-10°C MEDICAL REFRIGERATOR

Microprocessor controller, temperature ranging from +2°C to +8°C, can be set freely, controlling precision 1°C, display accuracy 0.1°C. Room temperature ranging from 0°C to 32°C

Combined with two exact sensors and auto defrost

Audible and visual alarm

Upright type, Exterior made from structural plate, interior made from ABS material,

Four units Caster are mounted under the bottom

Two-layer glass door, and inert gas inside. and lockable

Interior fluorescent lighting

Adjustable 4-5 units shelves made of quality steel wire

Forced air circulation system

Highly effective condenser and expansile evaporator to provide quick freezing

Optional: temperature printer or chart recorder



## SPECIFICATIONS

|                      |   |
|----------------------|---|
| Model                | REF11-300   |
| Capacity             | 300 L   |
| Temperature Range    | 2~8°C   |
| Refrigerant          | R134a, CFC free   |
| Refrigeration System | International famous compressor and Germany EBM fan motor           |
| Controller           | Microprocessor Control with LCD Displayer                           |
| Alarms               | High & low temperature alarm, Door open alarm, Sensor failure alarm |
| External Size        | 590x603x1265 mm   |
| Package Size         | 675x725x1625 mm   |
| Optional Accessories | Remote alarm system   |
| Weight               | 90/100 Kgs  |
| Power                | 215W  |
| Power Supply         | 220V,50/ 60HZ, 110V,50/60HZ   |



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)