

+4/+22 BLOOD TRANSPORT BOX REF21-090



+4/+22 BLOOD TRANSPORT BOX REF21-090

Designed as an essential piece of equipment in the immunohematology department and provides safe and convenient storage of whole blood, blood components, reagents and ensures freshness and integrity of blood and blood components. Microprocessor integrated system provides excellent temperature accuracy and reliable performance.

Used in Immunohematology Department, Blood Bank, Laboratory, Research, Medical, Pharmaceutical.

Also known as Laboratory Blood Bank Refrigerator.

REF21-090 +4/+22 BLOOD TRANSPORT BOX

4°C for storing blood

22°C for storing blood platelet

Audible and visual alarm

Four units of baskets inside as standard, which can load 50 blood bags of 450ml



SPECIFICATIONS

Model	REF21-090
Capacity	90 L
Temperature Range	4°C or 22°C
Ambient Temperature	-5°C or 43°C
Compressor	Danfoss compressor
Alarms	High & low temperature alarm
Insulation	50 mm
Material	Whole LLDPE
External Size	1000Lx565Wx862H mm
Package Size	1135Lx636Wx1069H mm
Optional Accessories	Hydraulic support bar, Remote alarm system
Weight	73/90 kg
Power Supply	DC12/24V/ AC100-240V 50/60Hz Automatic selection



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com