

# PORTABLE REFRIGERATOR REF31-240



# PORTABLE REFRIGERATOR REF31-240

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

## REF31-240 PORTABLE REFRIGERATOR



### SPECIFICATIONS

Model	REF31-240
External Size	330x235x275 mm
Internal Size	240x145x180 mm
Inner Capacity	6 L
Product Configuration	1pc refrigerator, 4 pcs Crystal storage ice box 400ml, 1pc digital thermometer
Insulation Material	PU
Gross Weight	2.8 kg
Packing Size	380x260x290 mm



### Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)