

PORTABLE REFRIGERATOR REF31-295



PORTABLE REFRIGERATOR REF31-295

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

REF31-295 PORTABLE REFRIGERATOR



SPECIFICATIONS

Model	REF31-295
External Size	335x225x205 mm
Internal Size	295x167x168 mm
Inner Capacity	8 L
Product Configuration	1pc refrigerator, 4 pcs Crystal storage ice box 400ml, 1pc digital thermometer
Insulation Material	EPS
Gross Weight	1.8 kg
Packing Size	380x260x290 mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com