

# PORTABLE REFRIGERATOR REF31-480



# PORTABLE REFRIGERATOR REF31-480

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

## REF31-480 PORTABLE REFRIGERATOR



### SPECIFICATIONS

Model	REF31-480
External Size	530x325x340 mm
Internal Size	480x270x265 mm
Inner Capacity	33 L
Product Configuration	1pc refrigerator, 12 pcs Crystal storage ice box 400ml 1pc digital thermometer
Insulation Material	PU
Gross Weight	5.5 kg
Packing Size	550x340x360 mm



#### Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)