

PORTABLE REFRIGERATOR REF32-420



PORTABLE REFRIGERATOR REF32-420

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

REF32-420 PORTABLE REFRIGERATOR

- 1.The storage container can keep content cooling up to 24 hours
- 2.Insulation material is PU CFC-free polyurethane
- 3.Made of environmental PP material,food grade, non-toxic
- 4.Storage container is durable for vaccine storage, blood products storage,dry ice storage,frozen specimen storag etc



SPECIFICATIONS

Model	REF32-420
External Size (WxDxH)	420x240x300 mm
Internal Size (WxDxH)	330x170x225 mm
Capacity	12 L
Temperature Range	2°C-8°C
Storage Time (Ambient 30°C)	24 h
Material	
External	PP
Internal	PU
Ice Vessels	7(CMB-400)
Max Loading	8 kg
Optional Accessory	LED display thermometer
Net Weight	2.5 kg
Gross Weight	6 kg
Package Size (WxDxH)	440x265x335 mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com