

# PORTABLE REFRIGERATOR REF32-570



# PORTABLE REFRIGERATOR REF32-570

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

## REF32-570 PORTABLE REFRIGERATOR

- 1.The storage container can keep content cooling up to 24 hours
- 2.Insulation material is PU CFC-free polyurethane
- 3.Made of environmental PP material,food grade, non-toxic
- 4.Storage container is durable for vaccine storage, blood products storage,dry ice storage,frozen specimen storag etc



## SPECIFICATIONS

|                             |                         |
|-----------------------------|-------------------------|
| Model                       | REF32-570               |
| External Size (WxDxH)       | 570x310x370 mm          |
| Internal Size (WxDxH)       | 450x225x300 mm          |
| Capacity                    | 33 L                    |
| Temperature Range           | 2°C-8°C                 |
| Storage Time (Ambient 30°C) | 20 h                    |
| Material                    |                         |
| External                    | PP                      |
| Internal                    | PU                      |
| Ice Vessels                 | 3(CMB-1200)             |
| Max Loading                 | 10 kg                   |
| Optional Accessory          | LED display thermometer |
| Net Weight                  | 4.5 kg                  |
| Gross Weight                | 9.5 kg                  |
| Package Size (WxDxH)        | 590x330x390 mm          |



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)