

PORTABLE REFRIGERATOR REF32-615



PORTABLE REFRIGERATOR REF32-615

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

REF32-615 PORTABLE REFRIGERATOR

- 1.The storage container can keep content cooling up to 24 hours
- 2.Insulation material is PU CFC-free polyurethane
- 3.Made of environmental PP material,food grade, non-toxic
- 4.Storage container is durable for vaccine storage, blood products storage,dry ice storage,frozen specimen storag etc



SPECIFICATIONS

| | |
|-----------------------------|-------------------------|
| Model | REF32-615 |
| External Size (WxDxH) | 615x395x450 mm |
| Internal Size (WxDxH) | 525x300x355 mm |
| Capacity | 54 L |
| Temperature Range | 2°C-8°C |
| Storage Time (Ambient 30°C) | 24 h |
| Material | |
| External | PP |
| Internal | PU |
| Ice Vessels | 4(CMB-1500) |
| Max Loading | 11 kg |
| Optional Accessory | LED display thermometer |
| Net Weight | 7.8 kg |
| Gross Weight | 17 kg |
| Package Size (WxDxH) | 800x550x510 mm |



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com