

# INCUBATOR SHAKER SHA29-0300



# INCUBATOR SHAKER SHA29-0300

In Incubator Shaker, the constant agitation with uniform amplitude provides the optimum oxygen transfer to the cultures, leading to better and consistent yield. A variety of platforms and different laboratory accessories permits configuration to fulfill your demands.

Used in Cell Cultures, Cell Aeration, Microbiology, Increasing Solubility Rates, Metabolism Studies, Bacterial Cultures, Bacteriology.

Also known as Laboratory Incubator Shaker, Incubator Shaker.

## SHA29-0300 INCUBATOR SHAKER

Microprocessor controller

Stainless steel interior and Electrostatic sprayed exterior

Flexible and convenient, wide variety of detachable accessories



## SPECIFICATIONS

Model	SHA29-0300
Speed Range	20-300 rpm
Standard Capacity	500mlx36
Temperature	RT+5-60°C
Temperature Accuracy	±0.1°C
Temperature Uniformity	±1°C
Timing Range	Optional
Motion	Orbital
Amplitude	26 mm
Tray	2
Platform Dimension	730x460 mm
Overall Dimension	980x820x1335 mm
Display	LCD
Weight	220 kg
Power	960 W
Power Supply	220V, 50/60Hz





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)