

# MICROPLATE SHAKER SHA32-1500



# MICROPLATE SHAKER SHA32-1500

Microplate Shaker is ideal for applications requiring precise speed control or pulse vortexing. Features true pulsing action for a more random motion and thorough mixing. It comes with interchangeable clamps for different micro plates which makes it universally fit.

Used in Stability, Dissolution Studies, Liquid Extractions, Protein Precipitation, Small Peptide Synthesis, Dilutions. Also known as Laboratory Orbital Shaker.

## SHA32-1500 MICROPLATE SHAKER

Variable speed from 200 rpm to 1,500 rpm

DC brushless motor for gentle and powerful vibration

Continuous or timing operation

Orbital of 3mm for efficient mixing in small sample vessels

Simple spring holders for easy and quick plate installation



## SPECIFICATIONS

Model	SHA32-1500
Speed Range	200~1500 rpm Standard & 200~1000 rpm Deep well
Speed Setting Resolution	10 rpm increment
Shaking Orbit	3 mm
Timing Range	1min~99h 59min
Overall Dimension	290x270x120 mm
Weight	6.5 kg
Power Supply	AC220V/AC120V, 2A



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)