

# MICROPLATE SHAKER SHA34-1500



# MICROPLATE SHAKER SHA34-1500

Microplate Shaker is ideal for applications requiring precise speed control or pulse vortexing. Features true pulsing action for a more random motion and thorough mixing. It comes with interchangeable clamps for different micro plates which makes it universally fit.

Used in Mixing, Incubation, Cultivation, Microbiology, Molecular biology, Biotechnology, Immunoassays, HTS applications, Colouration tests, Binding.

Also known as Laboratory Microplate Shaker, Microplate Shaker.

## SHA34-1500 MICROPLATE SHAKER

Easy to operate with one touch knob

DC brushless motor, long life and free of maintenance

Orbital of 3mm for efficient mixing in small sample vessels

Timing function



## SPECIFICATIONS

Model	SHA34-1500
Speed Range	200-1500 rpm
Capacity	4 microplates
Speed Increment	10 rpm
Shaking Orbit	3 mm
Timing Range	1 min ~ 99 h59 min
Overall Dimension	280×270×110 mm
Weight	6.5 kg
Power Supply	220 V/120 V, 2A



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)