

DECOLORISING SHAKER SHA62-25



DECOLORISING SHAKER SHA62-25

Rocking Shaker are designed to deliver salient performance, accuracy, and reproducible results for all your application needs. Lightweight and portable for easy transfer, outstanding uniform motion, high quality stain resistant platforms with non-slip rubber mats and low noise are the efficient and reliable features that makes it unique.

Used in Laboratory, Research, Mixing, Incubation.

Also known as Laboratory Rocking Shaker.

SHA62-25 DECOLORISING SHAKER

Rocking angle adjustable

Speed adjustable

Time setting



SPECIFICATIONS

Model	SHA62-25
Speed Range	0-25 rpm
Angle Range	±25°C
Inner Dimension	250×200 mm
Overall Dimension	420×200×200 mm
Weight	2.5 kg
Power Supply	AC 80-230 V 50/60 Hz



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com