

MINI MAGNETIC STIRRER STR11-250



MINI MAGNETIC STIRRER STR11-250

It is distinctly designed to meet the users application. Small, exquisite, convenient, portable, simple operation and long durability makes it an ideal choice.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

STR11-250 MINI MAGNETIC STIRRER

Distinct designed. Small and exquisite. Simple operation.

Long using life with no carbon brush motor.

Speed regulation stepless, large range of speed.

Strong magnetism to avoid stirrer jumping.

Convenient and portable.

Meet different kinds of experimental demands.



SPECIFICATIONS

Model	STR11-250
Capacity	Max. 250 ml
1.2	Motor rating output [W]
Stirrer Dimension	Φ8x40 mm
Platform Material	Plastic
Platform Dimension	Φ88 mm
Overall Dimension	125x110x70 mm
Weight	1.1 kg
Power	15 W
Power Supply	AC90~230V, 50/60 Hz



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com