

MAGNETIC STIRRER STR12-003



MAGNETIC STIRRER STR12-003

It is distinctly designed to meet the users application. Small, exquisite, convenient, portable, simple operation and long durability makes it an ideal choice.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

STR12-003 MAGNETIC STIRRER

Speed range of 0-1500 rpm

Maximum stirring quantity 3L

ABS casing with fire retardant material, resistance of weak acid and alkali



SPECIFICATIONS

Model	STR12-003
Capacity	3 L
Speed Range	0~1500 rpm
Temperature	5-40°C
Humidity	0.8
5	Motor rating output [W]
Work Plate Size	135 mm
Stirrer Dimension	50 mm
Platform Material	Plastic
Speed Display	Scale
Motor Type	DC motor
Protection Class	IP42
Overall Dimension	150Wx260Dx80H mm
Package Size (WxDxH)	220Wx295Dx180H mm
Weight	2 kg
Power	10 W
Power Supply	AC100/200V±10%, 50/60Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Dimension
LS55170	Stirrer bar		10mm x 6mm
LS55171	Stirrer bar		15mm x 8mm
LS55172	Stirrer bar		20mm x 8mm

LS55173	Stirrer bar		25mm x 8mm
LS55174	Stirrer bar		30mm x 6mm
LS55175	Stirrer bar		40mm x 8mm
LS55176	Stirrer bar		50mm x 8mm
LS55177	Stirrer bar remover	Length of 200mm	



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com