

MAGNETIC STIRRER STR13-001



MAGNETIC STIRRER STR13-001

It is distinctly designed to meet the users application. Small, exquisite, convenient, portable, simple operation and long durability makes it an ideal choice.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

STR13-001 MAGNETIC STIRRER

Maintenance free brushless DC motor guarantees long service

Rare earth magnets give powerful stirring

Stepless speed regulation

Ultra silence and vibration-free design.

PET top provides excellent chemical resistance

Shaft can be fitted with a pH sensor or sensors for other purposes



SPECIFICATIONS

| | |
|---------------------------------|-----------------------------|
| Model | STR13-001 |
| Capacity | 1.5 L |
| Stirring Positions | 1 |
| Speed Range | 300-2000 prm |
| Permissible Ambient Temperature | 5-40°C |
| Humidity | 80%RH |
| Work Plate Size | Φ120 mm |
| Stirrer Dimension | 40 mm |
| Platform Material | PET |
| Motor Type | Brushless DC motor |
| Protection Class | IP21 |
| Overall Dimension | 150Wx145Dx50H mm |
| Weight | 0.3 kg |
| Power | 7.2 W |
| Power Supply | 100-120/200-240 V, 50/60 Hz |



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com