

MAGNETIC STIRRER

STR14-008



MAGNETIC STIRRER STR14-008

It is distinctly designed to meet the users application. Small, exquisite, convenient, portable, simple operation and long durability makes it an ideal choice.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

STR14-008 MAGNETIC STIRRER

Speed range of 15-1500rpm

Automatic reverse rotation for better mixing result

High IP protection class (IP 65)

Set-up plate is made from chemically resistant materials

Slip-proof and stable



SPECIFICATIONS

Model	STR14-008
Capacity	0.8 L
Speed Range	15-1500 rpm
Temperature	5-40°C
Permissible Relative Humidity	0.8
5	Motor rating output [W]
Maximum Load Capacity	2 kg
Work Plate Size	90x90 mm
Platform Material	ABS
Protection Class	IP65
Overall Dimension	175Wx120Dx15H mm
Weight	0.3 kg
Power Supply	110-240V, 50/60 Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Dimension
LS55170	Stirrer bar		10mm x 6mm
LS55171	Stirrer bar		15mm x 8mm
LS55172	Stirrer bar		20mm x 8mm
LS55173	Stirrer bar		25mm x 8mm
LS55174	Stirrer bar		30mm x 6mm
LS55175	Stirrer bar		40mm x 8mm
LS55176	Stirrer bar		50mm x 8mm

LS55177	Stirrer bar remover	Length of 200mm	
---------	---------------------	-----------------	--



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com