

MAGNETIC STIRRER STR15-020



MAGNETIC STIRRER STR15-020

It is distinctly designed to meet the users application. Small, exquisite, convenient, portable, simple operation and long durability makes it an ideal choice.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

STR15-020 MAGNETIC STIRRER

Speed range of 100-1500 rpm

Max. stirring quantity 3L

LED display set and actual speeds

Microprocessor control technology guarantees speed control accuracy

ABS casing with fire retardant material, resistance of weak acid and alkali



SPECIFICATIONS

Model	STR15-020
Capacity	20 L
Speed Range	100~1500 rpm
Temperature	5-40°C
Humidity	0.8
5	Motor rating output [W]
Work Plate Size	135 mm
Stirrer Dimension	50 mm
Platform Material	Plastic
Speed Display	LED
Motor Type	DC motor
Protection Class	IP42
Overall Dimension	150Wx260Dx80H mm
Package Size (WxDxH)	220Wx295Dx180H mm
Weight	2 kg
Power	15 W
Power Supply	AC100/200V±10%, 50/60Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Dimension
LS55170	Stirrer bar		10mm x 6mm
LS55171	Stirrer bar		15mm x 8mm
LS55172	Stirrer bar		20mm x 8mm

LS55173	Stirrer bar		25mm x 8mm
LS55174	Stirrer bar		30mm x 6mm
LS55175	Stirrer bar		40mm x 8mm
LS55176	Stirrer bar		50mm x 8mm
LS55177	Stirrer bar remover	Length of 200mm	



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com