

MAGNETIC STIRRER STR17-010



MAGNETIC STIRRER STR17-010

It is distinctly designed to meet the users application. Small, exquisite, convenient, portable, simple operation and long durability makes it an ideal choice.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

STR17-010 MAGNETIC STIRRER

Speed range of 0-1500 rpm

Maximum stirring quantity 10 L

Glass ceramic plate, excellent chemical resistance and easy to clean



SPECIFICATIONS

Model	STR17-010
Capacity	10 L
Speed Range	0-1500 rpm

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Dimension
LS55170	Stirrer bar		10mm x 6mm
LS55171	Stirrer bar		15mm x 8mm
LS55172	Stirrer bar		20mm x 8mm
LS55173	Stirrer bar		25mm x 8mm
LS55174	Stirrer bar		30mm x 6mm
LS55175	Stirrer bar		40mm x 8mm
LS55176	Stirrer bar		50mm x 8mm
LS55177	Stirrer bar remover	Length of 200mm	
LS55195	Stirrer bar	Stirrer bar (65mm x 8mm), 1pc/pk	
LS55196	Stirrer bar	Stirrer bar (80mm x 13mm), 1pc/pk	



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com