

MULTISTATION MAGNETIC STIRRER STR32-08



MULTISTATION MAGNETIC STIRRER STR32-08

Multistation Magnetic Stirrer is designed for high-efficiency experiments and has stir bar configurations to choose. Laboratory space can be used effectively. Low speed but stable, high speed but powerful with low noise. Non-slip and anti-corrosion.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Multistation Magnetic Stirrer, Digital Magnetic Stirrer.

STR32-08 MULTISTATION MAGNETIC STIRRER

Digital display stirring speed at real-time, stir bar operate synchronously and provide uniform agitation reaction conditions.

Efficient Multistation magnetic stirrer configure 4/8/12 stirring positions to use lab space effective.

Low speed but stable, high speed but powerful. Low noise and maintenance free.

Stainless steel is coated by silicone resin film disk.

Non-slip and anti-corrosion.

DC brushless motor.



SPECIFICATIONS

Model	STR32-08
Capacity	400 ml
Stirring Positions	8
Speed Range	0~1100 rpm
Platform Material	Stainless steel with silicon resin film
Speed Display	Digital
Motor Type	Brushless DC motor
Overall Dimension	205Wx480Dx50H mm
Fuse	250V, 1A, Φ 5x20
Weight	4 kg
Power	20 W
Power Supply	AC100-230V, 50/60Hz



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com