

# MULTISTATION MAGNETIC STIRRER STR34-10



# MULTISTATION MAGNETIC STIRRER STR34-10

Multistation Magnetic Stirrer is designed for high-efficiency experiments and has stir bar configurations to choose. Laboratory space can be used effectively. Low speed but stable, high speed but powerful with low noise. Non-slip and anti-corrosion.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Multistation Magnetic Stirrer, Digital Magnetic Stirrer.

## STR34-10 MULTISTATION MAGNETIC STIRRER

Speed range of 0-1100 rpm

High performance magnetic stirrer with 10 stirring positions benefit laboratory space

Individually controlled stirring for consistency with various samples

Brushless DC motor, maintenance free and provides constant speed

Stainless steel plate with silicone film, anti-slip and anti-corrosion



## SPECIFICATIONS

|                                 |                                   |
|---------------------------------|-----------------------------------|
| Model                           | STR34-10                          |
| Capacity                        | 400 ml                            |
| Stirring Positions              | 10                                |
| Speed Range                     | 0~1100 rpm                        |
| Permissible Ambient Temperature | 5-40°C                            |
| Permissible Relative Humidity   | 0.8                               |
| Work Plate Size                 | 180x450 mm                        |
| Platform Material               | Stainless steel with silicon film |
| Speed Display                   | Scale                             |
| Motor Type                      | Brushless DC motor                |
| Overall Dimension               | 182Wx552Dx65H mm                  |
| Package Size (WxDxH)            | 250Wx710Dx165H mm                 |
| Weight                          | 7 kg                              |
| Power Supply                    | AC110/220V, 50/60Hz               |

## OPTIONAL ACCESSORIES

| Accessory Code | Name        | Description | Dimension  |
|----------------|-------------|-------------|------------|
| LS55170        | Stirrer bar |             | 10mm x 6mm |
| LS55171        | Stirrer bar |             | 15mm x 8mm |
| LS55172        | Stirrer bar |             | 20mm x 8mm |
| LS55173        | Stirrer bar |             | 25mm x 8mm |
| LS55174        | Stirrer bar |             | 30mm x 6mm |

|         |                     |                 |            |
|---------|---------------------|-----------------|------------|
| LS55175 | Stirrer bar         |                 | 40mm x 8mm |
| LS55177 | Stirrer bar remover | Length of 200mm |            |



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)