

# COAGULATION ANALYZER



# COAGULATION ANALYZER

Auto Coagulation Analyzer used to test clotting items, Chromogenic, Immunologic measuring methods  
Used in Clinical diagnostic, Laboratory, Pharmaceutical.  
Also known as coagulation analyzer.

## ANA13 AUTO COAGULATION ANALYZER



### SPECIFICATIONS

Model	ANA13-140	ANA13-160
Throughput	140T/H	160T/H
Analysis Method	Coagulation Method	Coagulation method, Immunospectrophotometry
Analysis Items	PT, APTT, TT, FIB	PT, APTT, TT, FIB, D-D, FDP
Sample Positions	5	6
Tilt reagent Positions	6	11 refrigerated reagent positions



ANA13-140



ANA13-160



### Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)