

COAGULATION ANALYZER



COAGULATION ANALYZER

Auto Coagulation Analyzer used to test clotting items, Chromogenic, Immunologic measuring methods
Used in Clinical diagnostic, Laboratory, Pharmaceutical.
Also known as coagulation analyzer.

ANA13 AUTO COAGULATION ANALYZER



SPECIFICATIONS

Model	ANA13-140	ANA13-160
Throughput	140T/H	160T/H
Analysis Method	Coagulation Method	Coagulation method, Immunoturbidimetry
Analysis Items	PT, APTT, TT, FIB	PT, APTT, TT, FIB, IDD, FDP
Sample Positions	5	6
Tilt reagent Positions	6	11 refrigerated reagent positions



ANA13-140



ANA13-160



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com