

PORTABLE REFRIGERATOR



PORTABLE REFRIGERATOR

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

REF31 PORTABLE REFRIGERATOR



SPECIFICATIONS

Model	REF31-240	REF31-295	REF31-325
External Size	330x235x275 mm	335x225x205 mm	425x262x305 mm
Internal Size	240x145x180 mm	295x167x168 mm	325x170x225 mm
Inner Capacity	6 L	8 L	12 L
Product Configuration	1pc refrigerator, 4 pcs Crystal storage ice box 400ml, 1pc digital thermometer		1pc refrigerator, 1pc tray 1pc partition 6 pcs Crystal storage ice box 400ml 1pc digital thermometer
Insulation Material	PU	EPS	PU
Gross Weight	2.8 kg	1.8 kg	4 kg
Packing Size	380x260x290 mm		450x330x330 mm

Model	REF31-363	REF31-370	REF31-480
External Size	425x262x305 mm	450x330x300 mm	530x325x340 mm
Internal Size	363x203x244 mm	370x235x240 mm	480x270x265 mm
Inner Capacity	15 L	20 L	33 L
Product Configuration	1pc refrigerator, 1pc strap 1pc tray 6 pcs Crystal storage ice box 400ml 1pc digital thermometer	1pc refrigerator, 1pc strap 1pc tray 9 pcs Crystal storage ice box 400ml 1pc digital thermometer 1pc partition	1pc refrigerator, 12 pcs Crystal storage ice box 400ml 1pc digital thermometer
Insulation Material	PU		
Gross Weight	4.5 kg	5 kg	5.5 kg
Packing Size	450x330x330 mm	490x350x330 mm	550x340x360 mm

Model	REF31-590	REF31-540	REF31-580
External Size	640x395x395 mm	620x485x375 mm	640x520x375 mm
Internal Size	590x335x320 mm	540x425x315 mm	580x465x320 mm

Inner Capacity	55 L	65 L	85 L
Product Configuration	1pc refrigerator, 12 pcs Crystal storage ice box 400ml 1pc digital thermometer	1pc refrigerator, 15 pcs Crystal storage ice box 400ml 1pc digital thermometer	
Insulation Material	PU		
Gross Weight	6.5 kg	8.2 kg	8.7 kg
Packing Size	650x410x410 mm	640x505x390 mm	660x540x395 mm



REF31-240



REF31-295



REF31-325



REF31-363



REF31-370



REF31-480



REF31-590



REF31-540



REF31-580

REF32 PORTABLE REFRIGERATOR

- 1.The storage container can keep content cooling up to 24 hours
- 2.Insulation material is PU CFC-free polyurethane
- 3.Made of environmental PP material,food grade, non-toxic
- 4.Storage container is durable for vaccine storage, blood products storage,dry ice storage,frozen specimen storag etc



SPECIFICATIONS

Model	REF32-340	REF32-420	REF32-440
External Size (WxDxH)	340x230x280 mm	420x240x300 mm	440x300x340 mm
Internal Size (WxDxH)	226x140x210 mm	330x170x225 mm	350x220x250 mm
Capacity	6 L	12 L	17 L
Temperature Range	2°C-8°C		
Storage Time (Ambient 30°C)	24 h		
Material			
External	PP		
Internal	PU		
Ice Vessels	6(CMB-400)	7(CMB-400)	2(CMB-1200)
Max Loading	6 kg	8 kg	15 kg
Optional Accessory	LED display thermometer		
Net Weight	1.9 kg	2.5 kg	3.2 kg
Gross Weight	5 kg	6 kg	6.5 kg
Package Size (WxDxH)	360x250x300 mm	440x265x335 mm	460x330x370 mm

Model	REF32-270	REF32-570	REF32-615
External Size (WxDxH)	440x270x310 mm	570x310x370 mm	615x395x450 mm
Internal Size (WxDxH)	375x215x250 mm	450x225x300 mm	525x300x355 mm
Capacity	20 L	33 L	54 L
Temperature Range	2°C-8°C		
Storage Time (Ambient 30°C)	18 h	20 h	24 h
Material			
External	PP		
Internal	PU		
Ice Vessels	8(CMB-400)	3(CMB-1200)	4(CMB-1500)
Max Loading	18 kg	10 kg	11 kg
Optional Accessory	LED display thermometer		
Net Weight	3.5 kg	4.5 kg	7.8 kg
Gross Weight	7.5 kg	9.5 kg	17 kg
Package Size (WxDxH)	460x290x330 mm	590x330x390 mm	800x550x510 mm



REF32-340



REF32-420



REF32-440



REF32-270



REF32-570



REF32-615



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com