

CHILDREN MEDICAL BED



CHILDREN MEDICAL BED

Children medical Bed built to the most current standards while also meeting the clinician's needs which allows easy care and disinfection. The bed is equipped with Trendelenburg and reverse Trendelenburg positions. Simple and modern design for safety and protection of the child.

Used in Medical unit, Nursing.

Also known as Children Hospital Bed, Medical bed, Hospital bed, Medical child bed, Pediatric bed, Medical pediatric bed, Hospital pediatric bed, Children medical bed, Children pediatric hospital bed.

HBD41-SS STAINLESS STEEL CHILDREN BED

- Stainless steel framework
- Stainless steel side rail, drop type
- Noiseless castors
- With mattress
- With debris basket



SPECIFICATIONS

Model	HBD41-SS
IV pole	Yes
Dimension	750x450x610-1030 mm

HBD42-SS1 STAINLESS STEEL CHILDREN BED

Stainless steel framework



SPECIFICATIONS

Model	HBD42-SS1
Dimension	1600x700x1000 mm

HBD44-03 MANUAL 3 FUNCTION CHILDREN MEDICAL BED

ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, pink, with patient name card holder, option: cartoon sticker

Full length powder coated side rail, drop type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

Central locking ABS foot pedal with casters dia. 5"

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes



SPECIFICATIONS

Model	HBD44-03
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	480~720 mm (without mattress)
IV pole	Yes
Dimension	1980Lx980Wx480~720 mm

HBD45-03 ELECTRIC 3 FUNCTION CHILDREN MEDICAL BED

ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, pink, with patient name card holder, option: cartoon sticker

Full length powder coated side rail, drop type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

Central locking S.S foot pedal with casters dia. 5"

With 3 motors from Taiwan

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes



SPECIFICATIONS

Model	HBD45-03
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	480~720 mm (without mattress)
IV pole	Yes
Dimension	1980Lx980Wx480~720 mm

HBD47-02 CHILDREN BED WITH TWO FUNCTION



SPECIFICATIONS

Model	HBD47-02
Load Bearing	150 kg
Dimension	1770x900 mm
Packing Dimension	wooden carton: 1880x980x1450 mm (assembled)
carton	
1880x1020x310 mm (detached)	Net Weight



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com