

DIALYSIS BED



DIALYSIS BED

Dialysis Bed is mainly used for setting purpose during dialysis treatment. This bed is best to be used when patients with end-stage renal disease spend many hours each week undergoing hemodialysis. We contribute toward enhanced care, through our innovative range of seating solutions. Our specialised dialysis bed, have been created with both patient comfort and nurse safety in mind.

Used in Dialysis treatment, Hemodialysis.

Also known as Hemodialysis bed, Medical bed, hospital bed, Hospital cot.

HBD72-02 MANUAL 2 FUNCTION DIALYSIS BED

ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder

2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish, 3 sections Bed framework made from mild steel ERW rectangular tubes

With casters dia. 5" with individual brake

With equipment shelf on the foot end

2 functions adjusted by hidden type hand crank

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessories: Shoe shelf, mosquito net rack, dinner table



SPECIFICATIONS

Model	HBD72-02
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2160Lx980Wx500H mm

HBD73-02 MANUAL 2 FUNCTION DIALYSIS BED

ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder

2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish, 3 sections Bed framework made from mild steel ERW rectangular tubes

Central locking ABS foot pedal with casters dia. 5"

With equipment shelf on the foot end

2 functions adjusted by hidden type hand crank

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessories: Shoe shelf, mosquito net rack, dinner table



SPECIFICATIONS

Model	HBD73-02
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2160Lx980Wx500H mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com