

MANUAL HOSPITAL BED



MANUAL HOSPITAL BED

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

HBD21-HM ABS FLAT BED



SPECIFICATIONS

Model	HBD21-HM
Dimension	2030x900x500 mm

HBD22-HM1 ABS BED SINGLE CRANK



SPECIFICATIONS

Model	HBD22-HM1
Dimension	2030x900x500 mm

HBD23 MANUAL HOSPITAL BED

ABS head & foot boards can be removable.

Aluminium alloy side rail, collapsed type with safe lock.

Bed sheet top made from perforated steel, epoxy powder coated and baked finish.

Bed framework made from mild steel ERW rectangular tubes.

Covered castors dia125mm with brake

With telescopic SS IV pole and prevision



SPECIFICATIONS

Model	HBD23-HM2	HBD23-HM3
Function	2	
Back-rest lifting angle	0°-80°±5°	
Foot-rest lifting angle	0°-45°±5°	
IV pole	Yes	
Dimension	2080x970x500 mm	



HBD23-HM2



HBD23-HM3

HBD24-MC MANUAL BED DOUBLE CRANK

ABS head and boards can be removable

ABS side rail, 4-piece, tuck away type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

Central locking foot pedal control with castors dia 150 mm

With shoes holder

With telescopic SS IV pole and previsions



SPECIFICATIONS

Model	HBD24-MC
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2080x970x500 mm

HBD25-MC1 MANUAL BED DOUBLE CRANK

ABS head & foot boards can be removable

Aluminium alloy side rail, collapsed type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

Covered castors dia125mm with brake

With drainage bag holders, with turning table and shoe holder

With telescopic SS IV pole and prevision



SPECIFICATIONS

Model	HBD25-MC1
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Drainage bag holders	Yes
IV pole	Yes
Dimension	2080x970x500 mm

HBD26-03 MANUAL BED 3 CRANK

ABS head and boards can be removable

ABS side rail, 4-piece, tuck away type With safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish, with soft connection.

Bed framework made from mild steel ERW rectangular tubes.

Central locking foot pedal control with castors dia 150 mm.

With drainage bag holders

With telescopic SS IV pole and previsions



SPECIFICATIONS

Model	HBD26-03
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	500-750 mm (without mattress)
Drainage bag holders	Yes
IV pole	Yes
Dimension	2180x1030x500 mm

HBD27-05 MANUAL BED FOUR CRANK 5 FUNCTION

ABS head & foot boards can be removable.

Aluminium alloy side rail, 6-staff, collapsed type with safe lock.

Bed sheet top made from perforated steel, epoxy powder coated and baked finish.

Bed framework made from mild steel ERW rectangular tubes.

Covered castors dia125mm with brake

With telescopic SS IV pole and prevision



SPECIFICATIONS

Model	HBD27-05
Function	5
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Trendelenburg	12°
Reverse trendelenburg	12°
Height adjustment	500-750 mm (without mattress)
IV pole	Yes
Dimension	2080x970x500 mm

HBD28-02 TWO CRANK MANUAL BED

Steel framework, top class ABS & PP plastic

Collapsible Aluminium alloy protective railing; Detachable PP head & foot board

4 pcs universal castor, dia. 125 mm, Individual brake

A shoe rack under the bed board



SPECIFICATIONS

Model	HBD28-02
Function	2
Back-rest lifting angle	0°-75°±5°
Foot-rest lifting angle	0°-45°±5°
Patient surface	1910Lx830W mm
Height adjustment	470 mm
IV pole	Yes
Dimension	2080Lx980Wx470 mm
Net Weight	65 kg

HBD29-03 3 CRANK MANUAL BED

Steel framework, top class ABS & PP plastic

PP head & foot board and bed board, detachable and easy to clean 6 columns collapsible Aluminium alloy protective railing

4 pcs universal castor, dia. 125 mm, Individual brake



SPECIFICATIONS

Model	HBD29-03
Function	3
Back-rest lifting angle	0°-65°±5°
Foot-rest lifting angle	0°-45°±5°
Patient surface	1930x830 mm
Height adjustment	750-575 mm

IV pole	Yes
Dimension	2160x960x530- 730 mm

HBD2A-04 4 CRANK MANUAL BED

Steel framework, top class ABS & PP plastic
 PP material head / foot board, detachable and easy-clean
 4 pcs drop-down PP protective railing, controlled by gas spring
 Central brake castors, dia. 150 mm. Brake pedal locate below the foot board



SPECIFICATIONS

Model	HBD2A-04
Function	5
Back-rest lifting angle	0°-60°±5°
Foot-rest lifting angle	0°-40°±5°
Trendelenburg	0-12°
Reverse trendelenburg	0-12°
Patient surface	1950x900 mm
Height adjustment	755-485 mm
IV pole	Yes
Dimension	2200Lx1050Wx485-755 mm
Net Weight	117 kg

HBD2B-FM FLAT MEDICAL BED

ABS head & foot boards can be removable. with bumper, 8kg/pair, color :blue, with patient name card holder
 2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock
 Bed sheet top made from perforated steel, epoxy powder coated and baked finish
 Bed framework made from mild steel ERW rectangular tubes
 4 pcs Individual lockable Taiwan Secure castors at 4 corners, diam 5"
 Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes
 Optional Accessories: Shoe shelf, mosquito net rack, dinner table



SPECIFICATIONS

Model	HBD2B-FM
Back-rest lifting angle	0°-75°
Foot-rest lifting angle	0°-35°
IV pole	Yes
Dimension	2160Lx980Wx500H mm

HBD2E-02 MANUAL 2 FUNCTION MEDICAL BED

ABS head & foot boards can be removable. with bumper, 8kg/pair, color option: brown, blue, with patient name card holder

2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock.

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

With casters dia. 5" with individual brake

Manual 2 functions adjusted by hidden type hand crank

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessories: Shoe shelf, mosquito net rack, dinner table



SPECIFICATIONS

Model	HBD2E-02
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2150Lx980Wx500H mm

HBD2G-03 MANUAL 3 FUNCTION MEDICAL BED

ABS head & foot boards can be removable. with bumper, 7kg/pair, color blue, with patient name card holder

2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

With casters dia. 5" with individual brake

Manual 3 functions adjusted by hidden type hand crank

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessories: Shoe shelf, mosquito net rack, dinner table



SPECIFICATIONS

Model	HBD2G-03
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	430~720 mm
IV pole	Yes
Dimension	2120Lx980Wx430~720 mm

HBD2H-02 MANUAL 2 FUNCTION MEDICAL BED

ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder

2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

Central locking ABS foot pedal with casters dia. 5"

Manual 2 functions adjusted by hidden type hand crank

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessories: Shoe shelf, mosquito net rack, dinner table



SPECIFICATIONS

Model	HBD2H-02
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2160Lx980Wx500 mm

HBD2J-02 MANUAL 2 FUNCTION MEDICAL BED

ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder

ABS side rail, 4-piece, tuck away type with safe lock. with angle display

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes.

Central locking both side foot pedal with casters dia. 5"

Manual 2 functions adjusted by hidden type hand crank

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes



SPECIFICATIONS

Model	HBD2J-02
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2220Lx1060Wx500 mm

HBD2K-03 MANUAL 3 FUNCTION MEDICAL BED

ABS head & foot boards can be removable.

with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder

ABS side rail, 4-piece, tuck away type with safe lock. with angle display

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes.

Central locking both side foot pedal with casters dia. 5".

Manual 3 functions adjusted by hidden type hand crank

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes



SPECIFICATIONS

Model	HBD2K-03
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	450~740 mm
IV pole	Yes
Dimension	2220Lx1060Wx450~740mm

HBD2M-HM SINGLE CRANK HOSPITAL BED



SPECIFICATIONS

Model	HBD2M-HM
Load Bearing	240 kg
Patient surface	2150x950 mm
Height adjustment	500 mm (from floor)

HBD2O-HM1 MANUAL HOSPITAL BED

Material: ABS

4 pcs of European style guardrail, upside down structure

Luxury central control castors, brake the bed with one stamp

With 4 bumpers



SPECIFICATIONS

Model	HBD2O-HM1
Back-rest lifting angle	0°-75°
Foot-rest lifting angle	0°-35°
Dimension	2280x1120x420-780 mm

HBD2P-HM2 MANUAL HOSPITAL BED WITH 2 FUNCTION



SPECIFICATIONS

Model	HBD2P-HM2
Back-rest lifting angle	0°-75°
Foot-rest lifting angle	0°-35°

Load Bearing	240 kg
Patient surface	1950x830 mm
Dimension	2200x1020x500 mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com